

Harmonic®

# Today's new standard for head and neck surgery

Precision + Fine dissection<sup>1</sup> + Improved efficiency<sup>2,3</sup>



HARMONIC FOCUS®+ Shears

**ETHICON**  
PART OF THE *Johnson & Johnson* FAMILY OF COMPANIES

Shaping  
the future  
of surgery

# Proven leader in thyroidectomy—now get the same **precision and speed in your head and neck procedures\***

The **HARMONIC FOCUS®+ Shears with Adaptive Tissue Technology** is the device surgeons need for fine dissection and sealing of vessels up to 5mm in all head and neck procedures. HARMONIC FOCUS®+ Shears enables safe and efficient surgery with less blood loss, delivering shorter procedures and better outcomes in the operating room.

**Speed:**  
Reduced total operative time by 31%<sup>4</sup>




**31%**  
(p<0.001)

**Outcomes:**  
Reduced intraoperative blood loss by 45ml<sup>5</sup>



**45**  
ml  
(p<0.001)

**Safety:**  
Same nerve function following use at 2mm from the sciatic nerve<sup>6</sup>




**EQUIVALENT**

**Cost-savings:**  
Reduced total operative costs by 10%<sup>7</sup>



**10%**  
(p=0.007)



Parotidectomy\*

Submandibular gland excision\*


Neck dissection\*

Parathyroidectomy\*

Thyroidectomy\*

\* Approved ENT indication, see IFU

# HARMONIC FOCUS<sup>®</sup>+ Shears— designed for **efficiency, fine dissection and precision**



## ORDERING INFORMATION

### HARMONIC FOCUS<sup>®</sup>+ Shears with Adaptive Tissue Technology

Product code	Description	Quantity per sales unit
HAR9F	HARMONIC FOCUS <sup>®</sup> + Shears with Adaptive Tissue Technology	6

To learn more about Ethicon solutions for head and neck, contact your sales representative or visit [ethicon.com/focus](http://ethicon.com/focus).

# Providing a comprehensive solution for head and neck procedures



**HARMONIC FOCUS®+ Shears** delivers fine dissection and sealing of vessels up to 5mm in all head and neck procedures

## Addressing surgical bleeding with adjunctive hemostats

**SURGICEL Family of Absorbable Hemostats**—50-plus years of the proven safety and efficacy surgeons trust.<sup>12,13</sup>

SURGICEL® FIBRILLAR™ lightweight layered structure is customizable and “melts” into bleeding tissue; a versatile option for use at multiple sites<sup>14</sup>



## Secure closure to address a known risk factor associated with infection

**STRATAFIX™ Symmetric PDS™ Plus Knotless Tissue Control Devices**—deliver more consistency, more security and more efficiency than traditional sutures<sup>15,16,17,18,19,20,21</sup>

**DERMABOND® PRINEO® Skin Closure System**—uncompromised strength and protection for excellent wound closure<sup>22</sup>



**1** (C260) **2** Based on a meta-analysis of HARMONIC FOCUS® (HF) versus clamp, cut and tie, where HF reduced operative time (p<0.001), intra-operative blood loss (p<0.001), length of stay (p<0.005), drainage volume (p<0.01). Cheng et al. A systematic review and meta-analysis of Harmonic Focus in thyroidectomy compared to conventional techniques. *Thyroid Research* (2015) 8:15 (C1962) **3** The health technology method was applied in a case study of 440 patients undergoing thyroidectomy in Terni, Italy. The use of HARMONIC FOCUS® resulted in reducing overall procedure time from 143.33 minutes to 113.7 minutes (20.67%) and reducing overall hospital cost from €3,055 to €2,768 (9.39%). Lucchini R, et al. Health technology assessment and thyroid surgery. *Il Giornale di Chirurgia* (July/August 2013) 34:198-201. (C1529) **4** Based on a meta-analysis of HARMONIC FOCUS® (HF) versus clamp, cut and tie, where HF reduced operative time (p<0.001). Cheng et al. A systematic review and meta-analysis of Harmonic Focus in thyroidectomy compared to conventional techniques. *Thyroid Research* (2015) 8:15. **5** Based on a meta-analysis of HARMONIC FOCUS® (HF) versus clamp, cut and tie, where HF reduced intra-operative blood loss (p<0.001). Cheng et al. A systematic review and meta-analysis of Harmonic Focus in thyroidectomy compared to conventional techniques. *Thyroid Research* (2015) 8:15. **6** In a preclinical model, there was no difference in nerve function following the use of HARMONIC FOCUS® shears at 2mm from the sciatic nerve versus cold scissors. **7** Cheng H et al. Hospital costs associated with thyroidectomy performed with a Harmonic device compared to conventional techniques: a systematic review and meta-analysis. *J Med Econ*. 2016 Apr 5:1-9. [Epub ahead of print] **8** (C260) **9** Based on preclinical testing of lymphatic vessels up to and including 1mm in diameter. **10** (C301) **11** (C1568) **12** A Compendium of Scientific Literature. Evidence Supporting the Efficacy and Safety of the Surgicel Family of Absorbable Hemostats **13** Hong Y, Loughlin K. The use of hemostatic agents and sealants in urology. *J Urol*. 2006;176:2367-2374. **14** Harrington, Claudia G. Hemostatic effectiveness of SURGICEL®, SURGICEL NU-KNIT®, and SURGICEL® FIBRILLAR™ in Swine Spleen Model. Ethicon, Inc. Data on file. 2005 Aug 1. **15** Eickmann T, Quane E. Total knee arthroplasty closure with barbed sutures. *J Knee Surg*. 2010;23(3):163-167. **16** Einarsson JJ, Chavan NR, Suzuki Y, Jonsdottir G, Vellinga TT, Greenberg JA. Use of bidirectional barbed suture in laparoscopic myomectomy: evaluation of perioperative outcomes, safety, and efficacy. *J Minim Invasive Gynecol*. 2011;18(1):92-95. **17** Levine BR, Ting N, Della Valle CJ. Use of a barbed suture in the closure of hip and knee arthroplasty wounds. *Orthopedics*. 2011;34(9):e473-e475. doi: 10.3928/01477447-20110714-35. **18** Moran ME, Marsh C, Perrotti M. Bidirectional-barbed sutured knotless running anastomosis v classic Van Velthoven suturing in a model system. *J Endourol*. 2007;21(10):1175-1178. **19** Rodeheaver GT, Pineros-Fernandez A, Salopek LS, et al. Barbed sutures for wound closure: in vivo wound security, tissue compatibility and cosmesis measurements. In: Transactions from the 30th Annual Meeting of the Society for Biomaterials; Mount Laurel, NJ; p. 232. **20** Vakili JJ, O'Reilly MP, Sutter EG, Mears SC, Belkoff SM, Khanuja HS. Knee arthroscopy repair with a continuous barbed suture: a biomechanical study. *J Arthroplasty*. 2011;26(5):710-713. **21** Warner JP, Gutowski KA. Abdominoplasty with progressive tension closure using a barbed suture technique. *Aesthet Surg J*. 2009;29(3):221-225. **22** Data on file, Ethicon, Inc.; DERMABOND PRINEO Skin Closure System Matrix. PRI-04413. 2013.