


PROXIMATE® Reloadable Linear Staplers (TL, TX)

Parallel jaw closure designed for consistent staple formation



Parallel jaw closure helps ensure that tissue is compressed evenly within the jaws of the instrument for consistent staple formation

- Two staggered rows of titanium alloy staples in 30mm, 60mm, 90mm sizes
- Three staggered rows of titanium alloy staples in 30mm PROXIMATE TX30V

Ethicon reloadable linear staplers allow the surgeon to choose the staple height to better accommodate various tissue thicknesses

PROXIMATE TL

- Dial allows controlled tissue compression and adjustable staple height for tissue management
- Maximum of 4 firings

PROXIMATE TX

- Three different reloads accommodate a variety of tissue thicknesses
- Maximum of 8 firings

PROXIMATE TL	
Feature	Benefit
Adjustable Staple Height	Surgeon controls staple height for tissue management
Heavy Wire option	Ability to staple thick tissue

PROXIMATE TX	
Feature	Benefit
Intermediate Locking Position	Facilitates proper placement of stapler on tissue
One handed operation (close, fire, release)	Easy to use

Applicable specialties

- Alimentary tract, Thoracic

Ordering Information

PROXIMATE® Reloadable Linear Stapler (TL)

Device Code	Staple Rows	Staple Line Length	Number of Staples	Qty
TLV30	2	30mm	15	3
TL30	2	30mm	11	3
TLH30	2	30mm	11	3
TL60	2	60mm	21	3
TLH60	2	60mm	21	3
TL90	2	90mm	33	3
TLH90	2	90mm	33	3

Reload Code	Closed Staple Height (mm)	Qty
TRV30	Secures Tissue Thickness of 1.0mm	12
TR30	Secures Tissue Thickness of 1.0 - 2.5mm	12
TRH30	Secures Tissue Thickness of 1.0 - 2.5mm	12
TR60	Secures Tissue Thickness of 1.0 - 2.5mm	12
TRH60	Secures Tissue Thickness of 1.5 - 2.5mm	12
TR90	Secures Tissue Thickness of 1.0 - 2.5mm	12
TRH90	Secures Tissue Thickness of 1.5 - 2.5mm	12

PROXIMATE® Reloadable Linear Stapler (TX)

Device Code	Staple Rows	Staple Line Length	Number of Staples	Qty
TX30V	3	30mm	23	3
TX30B	2	30mm	11	3
TX30G	2	30mm	11	3
TX60B	2	60mm	21	3
TX60G	2	60mm	21	3

Reload Code	Closed Staple Height (mm)	Qty
XR30V	1.0mm	12
XR30B	1.5mm	12
XR30G	2.0mm	12
XR60B	1.5mm	12
XR60G	2.0mm	12

How to Order

All purchase orders are made to Johnson & Johnson Health Care Systems, Inc. (JJHCS) or your authorized distributor.

Electronic Ordering Options

(Note: Placing orders electronically avoids minimum order fees for your hospital.)

Johnson & Johnson Gateway:

www.jnjgateway.com/commerce

For questions about your order, visit the website or call 1-866-JNJ-GATE (565-4283).

Global Healthcare Exchange (GHX):

www.ghx.com

For questions about your order, visit the website or call 1-800-YOUR-GHX (968-7449).

Electronic Data Interchange (EDI):

JJHCS EDI Help Line: 1-800-262-2888

Nonelectronic/Manual Ordering Options

Call JJHCS at 1-800-255-2500 between 8:30 a.m. and 9:00 p.m. Eastern Standard Time or fax your order to 1-732-562-2212.

Customer support

For product use assistance, clinical guidelines, service and repair, emergency assistance, copy of a 510(k) clearance letter or complaints, please contact our Customer Support Center at customersupport@eesus.jnj.com or by calling **1-877-ETHICON** (384-4266). Our Support Center is staffed 24 hours a day, 7 days a week by qualified nurses to answer your product-related questions.

For complete product details, see Instructions for Use.